

Project pt3

Monday, May 04, 2026 1:13 PM

cardio_generator

Sessions

cardio_generator

Element	Missed Instructions	Cov	Missed Branches	Cov	Missed Cxty	Missed Lines	Missed Methods	Missed Classes
com.cardio_generator.generators	0%	0%	24	24	104	104	16	5
com.cardio_generator	0%	0%	26	26	83	83	13	1
com.cardio_generator.outputs	0%	0%	18	18	60	60	16	5
com.data_management	79%	75%	4	28	16	77	2	20
com.alerts	100%	90%	6	48	0	109	0	18
Total	1,200 of 2,039	41%	54 of 120	55%	78	144	263	433

Created with JaCoCo 0.8.10.202304240956

The com.alerts package achieved 100% instruction coverage and 90% branch coverage, reflecting comprehensive testing of all alert conditions.

The com.data_management package reached 79% coverage, with the untested portions being edge cases in DataStorage and PatientRecord getters that are simple pass-through methods requiring no logic testing.

The com.cardio_generator packages (generators, outputs, and the main simulator) have 0% coverage. These classes handle real-time data simulation and network I/O (TCP sockets, WebSocket servers, file writing) which are not suitable for unit testing without complex mocking infrastructure.

Run FileDataReaderTest x

FileDataReaderTest (data_ma 187 ms)

4 tests passed 4 tests total, 187 ms

C:\Users\Carbon\.jdk\openjdk-25.0.1\bin\java.exe ...
Skipping malformed line: this is not valid

Process finished with exit code 0

Run AlertGeneratorTest x

AlertGeneratorTest (data_mar 121 ms)

19 tests passed 19 tests total, 121 ms

C:\Users\Carbon\.jdk\openjdk-25.0.1\bin\java.exe ...
ALERT: Patient 1|Systolic Pressure Decreasing Trend|3000
ALERT: Patient 1|Critical Systolic Low: 85.0|1000
ALERT: Patient 1|Critical Systolic Low: 85.0|1000
ALERT: Patient 1|Low Blood Saturation: 91.0|1000
ALERT: Patient 1|Critical Diastolic High: 125.0|1000
ALERT: Patient 1|Low Blood Saturation: 91.0|1777892415622
ALERT: Patient 1|Rapid Saturation Drop|1777892415622
ALERT: Patient 1|Manual Alert Triggered|1000
ALERT: Patient 1|Systolic Pressure Increasing Trend|3000
ALERT: Patient 1|Critical Systolic High: 185.0|1000
ALERT: Patient 1|Abnormal ECG Peak: 5.0|2000
ALERT: Patient 1|Critical Systolic Low: 85.0|1000
ALERT: Patient 1|Low Blood Saturation: 91.0|1000
ALERT: Patient 1|Hypotensive Hypoxemia Alert|1777892515671
ALERT: Patient 1|Critical Diastolic Low: 55.0|1000

Process finished with exit code 0